



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Bromocresol Green-Methyl Red
Catalog:	IS10786
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Bromocresol Green, Sodium Salt	62625-32-5	ACS
Methyl Red, Sodium Salt	845-10-3	ACS

Specification:

Test	Specification
Appearance	Dark colored liquid
Suitability for Use	Greenish-Blue (pH 5.4) – Light Pink (pH 4.5), with Gray warning color

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

