



Product Specification

Product Information:

Products:	Murexide 0.25% in Ammonium Chloride
Catalog:	IS10793
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Ammonium Chloride	12125-02-9	ACS
Murexide	3051-09-0	ACS

Specification:

Test	Specification
Appearance	Rose-pink solid
Suitability for Use	Pink (Excess Calcium) - Purple(Excess EDTA)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

