



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	P.A.N. 0.1% in Methanol
Catalog:	IS10798
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Methyl Alcohol	67-56-1	ACS
1,2-Pyridylazo-2-naphthol	85-85-8	Indicator

Specification:

Test	Specification
Appearance	Amber liquid
Suitability for Use	Deep Red (Uncomplexed Copper Ions) - Yellow (Excess EDTA)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

