



Product Specification

Product Information:

Products:	Phenol Red 0.02%
Catalog:	IS10800
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Phenol Red, Sodium Salt	34487-61-1	ACS

Specification:

Test	Specification
Appearance	Red-orange liquid
pH at 25°C	7.15-7.25
Suitability for Use	Yellow (pH <6.8)-Pinkish-Orange (pH 7.0) - Pink(pH >8.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

