



Product Specification

Product Information:

Products:	Phenolphthalein 0.5% in 50% Isopropanol
Catalog:	IS10805
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Isopropyl Alcohol	67-63-0	ACS/USP
Phenolphthalein	77-09-8	ACS
Sodium Hydroxide	1310-73-2	Reagent

Specification:

Test	Specification
Appearance	Colorless to slightly pink liquid
Suitability for Use	Colorless (pH 8.0) - Pink (pH 10.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

