



Product Specification

Product Information:

Products:	Phenolphthalein 1% in Isopropanol
Catalog:	IS10807
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Isopropyl Alcohol	67-63-0	ACS
Phenolphthalein	77-09-8	ACS
Water	7732-18-5	ACS/ASTM/USP/EP
Sodium Hydroxide	1310-73-2	Reagent

Specification:

Test	Specification
Appearance	Slightly pinkish-gold liquid
Suitability for Use	Colorless (pH < 8.0) -Pink (pH >10.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

