



# Product Specification

## Product Information:

<b>Products:</b>	Starch 0.25% with Salicylic Acid
<b>Catalog:</b>	IS10809
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Starch, soluble	9005-84-9	ACS
Salicylic Acid	69-72-7	ACS

## Specification:

Test	Specification
Appearance	White translucent liquid
Suitability for Use	Colorless (Iodine absent) – Blue (Iodine present)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

