



Tri Stains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

# Product Specification

## Product Information:

|                  |                    |
|------------------|--------------------|
| <b>Products:</b> | Thymol Blue 0.04%  |
| <b>Catalog:</b>  | IS10812            |
| <b>Grade:</b>    | Indicator solution |

## Reagent Used:

| Reagents                 | CAS #      | Grade           |
|--------------------------|------------|-----------------|
| Water                    | 7732-18-5  | ACS/ASTM/USP/EP |
| Thymol Blue, Sodium Salt | 62625-21-2 | ACS             |

## Specification:

| Test                | Specification  |
|---------------------|--|
| Appearance          | Dark green liquid  |
| Suitability for Use | Red (pH < 1.2) - Yellow (pH 2.8 to pH 8.0) - Blue (pH > 9.2) |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

