



Product Specification

Product Information:

Products:	Xylenol Orange 0.1%
Catalog:	IS10816
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Xylenol Orange Sodium Salt	3618-43-7	ACS

Specification:

Test	Specification
Appearance	Reddish Orange to Dark Purple
Suitability for Use	Purple (Excess Metal) – Yellow (Excess EDTA)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

