



Product Specification

Product Information:

Products:	Crystal Violet Indicator, 1.0 % in Glacial Acetic Acid
Catalog:	IS10829
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Acetic Acid	64-19-7	ACS
Gentian Violet	548-62-9	ACS

Specification:

Test	Specification
Appearance	Dark blue to violet/purple liquid
Suitability for Use	Purple (Basic) - Blue-Green (Acidic)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

