



# Product Specification

## Product Information:

<b>Products:</b>	Congo Red Indicator 0.1%
<b>Catalog:</b>	IS10833
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Congo Red	573-58-0	Indicator

## Specification:

Test	Specification
Appearance	Red liquid
Suitability for Use	Red (pH > 5.0) - Blue (pH < 3.0)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

