



# Product Specification

## Product Information:

<b>Products:</b>	Diphenyl carbazone-Bromophenol Blue
<b>Catalog:</b>	IS10834
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Ethyl Alcohol	64-17-5	Reagent
Methyl Alcohol	67-56-1	Reagent
s-Diphenyl carbazone	538-62-5	ACS
Diphenyl carbazide	140-22-7	ACS
Bromophenol Blue, Sodium Salt	62625-28-9	ACS

## Specification:

Test	Specification
Appearance	Deep red liquid
Suitability for Use	Yellow (excess Chloride & proper pH) - Purple (excess Mercury)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

