



Product Specification

Product Information:

Products:	Methyl Orange 0.05%
Catalog:	IS10836
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Methyl Orange, Sodium Salt	547-58-0	ACS

Specification:

Test	Specification
Appearance	Orange liquid
Suitability for Use	Yellow (pH >4.4) - Pinkish-Orange (pH 3.2) – Pink (pH < 3.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

