



Product Specification

Product Information:

Products:	Methyl Red 0.02% In Methanol
Catalog:	IS10839
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Methyl Alcohol	67-56-1	ACS
Methyl Red	493-52-7	ACS

Specification:

Test	Specification
Appearance	Red liquid
Suitability for Use	Pink (pH < 4.2) - Yellow (pH > 6.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

