



Product Specification

Product Information:

Products:	Methyl Red-Methylene Blue
Catalog:	IS10840
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Ethyl Alcohol	64-17-5	Reagent
Methyl Alcohol	67-56-1	Reagent
Methyl Red, Sodium Salt	845-10-3	ACS
Methylene Blue Trihydrate	7220-79-3	Certified

Specification:

Test	Specification
Appearance	Dark green liquid
Suitability for Use	Purple (pH < 4.2) - Green (pH > 6.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

