



Tri Stains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

# Product Specification

## Product Information:

<b>Products:</b>	Pan Indicator 0.3% In Etoh
<b>Catalog:</b>	IS10842
<b>Grade:</b>	Indicator solution

## Reagent Used:

Reagents	CAS #	Grade
Ethyl Alcohol	64-17-5	ACS
1,2-Pyridylazo-2-naphthol	85-85-8	ACS

## Specification:

Test	Specification
Appearance	Amber liquid
Suitability for Use	Deep Red (Uncomplexed Copper Ions) - Yellow (Excess EDTA)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

