

Product Specification

Product Information:

Products:	Starch Indicator 0.5% (Preserved with Chloroform)
Catalog:	IS10847
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Starch, soluble	9005-84-9	ACS
Chloroform	67-66-3	

Specification:

Test	Specification
Appearance	Translucent liquid
Suitability for Use	Colorless (Iodine absent) – Blue (Iodine present)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

