



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product Information:

Products:	Starch Indicator, 0.2%
Catalog:	IS10851
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Starch, soluble	9005-84-9	ACS
Salicylic Acid	69-72-7	ACS

Specification:

Test	Specification
Appearance	White translucent liquid
Suitability for Use	Colorless (Iodine absent) – Blue (Iodine present)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

