



Product Specification

Product Information:

Products:	Thymol Blue 0.4%
Catalog:	IS10854
Grade:	Indicator solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/ASTM/USP/EP
Thymol Blue, Sodium Salt	62625-21-2	ACS

Specification:

Test	Specification
Appearance	Dark green liquid
Suitability for Use	Red (pH < 1.2) - Yellow (pH 2.8 to pH 8.0) - Blue (pH > 9.2)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

