



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Fuchsin Basic Certified
Catalog# : S10092
CAS : 632-99-5
CI number : 42510
Formula : $C_{20}H_{20}N_3Cl$
MW : 337.86

TEST	SPECIFICATION
Appearance	Dark Green Crystalline Powder
Identification	To pass tests
Loss on drying at 110°C/1hour	$\leq 10.0\%$
Dye content	$\geq 88.0\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

