



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : New Coccine
Catalog# : S10127
CAS : 2611-82-7
CI number : -
Formula : C₂₀H₁₁N₂Na₃O₁₀S₃
MW : 604.47

TEST	SPECIFICATION
Appearance	Dark red powder
Insoluble	0.072% (in water)
Solubility	Clear, red to dark red solution (1 mg/mL in water)
Dye Content	84.21%
Heavy Metals	< 40 ppm
Arsenate	< 3 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

